

Product Datasheet

Human Chemokine C-X3-C-Motif Receptor 1 (CX3CR1) ELISA Kit (orb1207117)

Description	Human Chemokine C-X3-C-Motif Receptor 1 (CX3CR1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P49238
Sensitivity	0.055ng/mL
Entrez	1524
Expiration Date	Please enquire.