

Product Datasheet

Human Cholinergic Receptor, Nicotinic, Beta 2 (CHRNb2) ELISA Kit (orb1207131)

Description	Human Cholinergic Receptor, Nicotinic, Beta 2 (CHRNb2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P17787
Sensitivity	0.054ng/mL
Entrez	1141
Expiration Date	Please enquire.