

Product Datasheet

Human CUB Domain Containing Protein 1 (CDCP1) ELISA Kit (orb1207146)

Description	Human CUB Domain Containing Protein 1 (CDCP1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9H5V8
Sensitivity	0.127ng/mL
Entrez	64866
Expiration Date	Please enquire.