

Product Datasheet

Human Chemokine C-C-Motif Receptor 7 (CCR7) ELISA Kit (orb1207159)

Description	Human Chemokine C-C-Motif Receptor 7 (CCR7) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P32248
Sensitivity	0.059ng/mL
Entrez	1236
Expiration Date	Please enquire.