

Product Datasheet

Human Calpain 1, Large Subunit (CAPN1) ELISA Kit (orb1207167)

Description	Human Calpain 1, Large Subunit (CAPN1) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P07384
Sensitivity	0.58ng/mL
Entrez	823
Expiration Date	Please enquire.