

Product Datasheet

Human Calpain 10 (CAPN10) ELISA Kit (orb1207168)

Description	Human Calpain 10 (CAPN10) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9HC96
Sensitivity	0.047ng/mL
Entrez	11132
Expiration Date	Please enquire.