

Product Datasheet

Human Bone Marrow Stromal Cell Antigen 1 (BST1) ELISA Kit (orb1207185)

Description	Human Bone Marrow Stromal Cell Antigen 1 (BST1) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q10588
Sensitivity	33pg/mL
Entrez	683
Expiration Date	Please enquire.