

Product Datasheet

Human Beta-1,4-Galactosyltransferase 1 (b4GALT1) ELISA Kit (orb1207194)

Description	Human Beta-1,4-Galactosyltransferase 1 (b4GALT1) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P15291
Sensitivity	0.57ng/mL
Entrez	2683
Expiration Date	Please enquire.