

Product Datasheet

Human Arachidonate-5-Lipoxygenase (ALOX5) ELISA Kit (orb1207229)

Description	Human Arachidonate-5-Lipoxygenase (ALOX5) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P09917
Sensitivity	0.057ng/mL
Entrez	240
Expiration Date	Please enquire.