

Product Datasheet

Human Arachidonate-12-Lipoxygenase (ALOX12) ELISA Kit (orb1207230)

Description	Human Arachidonate-12-Lipoxygenase (ALOX12) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P18054
Sensitivity	0.116ng/mL
Entrez	239
Expiration Date	Please enquire.