

Product Datasheet

Human Adipose Differentiation Related Protein (ADRP) ELISA Kit (orb1207239)

Description	Human Adipose Differentiation Related Protein (ADRP) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q99541
Sensitivity	0.057ng/mL
Entrez	123
Expiration Date	Please enquire.