

Product Datasheet

Human Adenylate Cyclase Activating Polypeptide 1, Pituitary (ADCYAP1) ELISA Kit (orb1207241)

Description	Human Adenylate Cyclase Activating Polypeptide 1, Pituitary (ADCYAP1) ELISA Kit
Reactivity	Human
Range	7.81-500pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P18509
Sensitivity	2.3pg/mL
Entrez	116
Expiration Date	Please enquire.