

Product Datasheet

Human Adenylate Cyclase 1, Brain (ADCY1) ELISA Kit (orb1207244)

Description	Human Adenylate Cyclase 1, Brain (ADCY1) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q08828
Sensitivity	0.28ng/mL
Entrez	107
Expiration Date	Please enquire.