

Product Datasheet

Human Actin Gamma 2, Smooth Muscle (ACTg2) ELISA Kit (orb1207251)

Description	Human Actin Gamma 2, Smooth Muscle (ACTg2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P63267
Sensitivity	0.053ng/mL
Entrez	72
Expiration Date	Please enquire.