

Product Datasheet

Human ATP Binding Cassette Transporter A1 (ABCA1) ELISA Kit (orb1207258)

Description	Human ATP Binding Cassette Transporter A1 (ABCA1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O95477
Sensitivity	0.053ng/mL
Entrez	19
Expiration Date	Please enquire.