

Product Datasheet

Rat Platelet Derived Growth Factor Subunit B (PDGFB) ELISA Kit (orb1207266)

Description	Rat Platelet Derived Growth Factor Subunit B (PDGFB) ELISA Kit
Reactivity	Rat
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q05028
Sensitivity	12.5pg/mL
Entrez	24628
Expiration Date	Please enquire.