

Product Datasheet

Rat Interleukin 22 Receptor Alpha 2 (IL22Ra2) ELISA Kit (orb1207272)

Description	Rat Interleukin 22 Receptor Alpha 2 (IL22Ra2) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q7TNI4
Sensitivity	0.105ng/mL
Entrez	444986
Expiration Date	Please enquire.