

Product Datasheet

Rat Interleukin 20 (IL20) ELISA Kit (orb1207273)

Description	Rat Interleukin 20 (IL20) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	D4A2T8
Sensitivity	5.7pg/mL
Entrez	681112
Expiration Date	Please enquire.