

Product Datasheet

Rat Fibroblast Growth Factor 3 (FGF3) ELISA Kit (orb1207277)

Description	Rat Fibroblast Growth Factor 3 (FGF3) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q8R5L9
Sensitivity	5.8pg/mL
Entrez	170633
Expiration Date	Please enquire.