

Product Datasheet

Rat Endothelin 2 (EDN2) ELISA Kit (orb1207284)

Description	Rat Endothelin 2 (EDN2) ELISA Kit
Reactivity	Rat
Range	7.81-500pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P23943
Sensitivity	2.96pg/mL
Entrez	24324
Expiration Date	Please enquire.