

Product Datasheet

Porcine Vascular Endothelial Growth Factor A (VEGFA) ELISA Kit (orb1207288)

Description	Porcine Vascular Endothelial Growth Factor A (VEGFA) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P49151
Sensitivity	6.3pg/mL
Entrez	397157
Expiration Date	Please enquire.