

Product Datasheet

Bovine Transforming Growth Factor Beta 1 (TGFB1) ELISA Kit (orb1207296)

Description	Bovine Transforming Growth Factor Beta 1 (TGFB1) ELISA Kit
Reactivity	Bovine
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P18341
Sensitivity	10.7pg/mL
Entrez	282089
Expiration Date	Please enquire.