

Product Datasheet

Bovine Interleukin 6 (IL6) ELISA Kit (orb1207302)

Description	Bovine Interleukin 6 (IL6) ELISA Kit
Reactivity	Bovine
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P26892
Sensitivity	31pg/mL
Entrez	280826
Expiration Date	Please enquire.