

Product Datasheet

Bovine Interleukin 5 (IL5) ELISA Kit (orb1207304)

Description	Bovine Interleukin 5 (IL5) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P52173
Sensitivity	6.2pg/mL
Entrez	280825
Expiration Date	Please enquire.