

## Product Datasheet

### Bovine Interleukin 18 (IL18) ELISA Kit (orb1207316)

|                          |                                                                                                            |
|--------------------------|------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Bovine Interleukin 18 (IL18) ELISA Kit                                                                     |
| <b>Reactivity</b>        | Bovine                                                                                                     |
| <b>Range</b>             | 15.6-1000pg/mL                                                                                             |
| <b>Note</b>              | For research use only                                                                                      |
| <b>Application notes</b> | sandwich                                                                                                   |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates,<br>cell lysates, cell culture supernates or<br>other biological fluids. |
| <b>Uniprot ID</b>        | <b>Q9TU73</b>                                                                                              |
| <b>Sensitivity</b>       | 6.1pg/mL                                                                                                   |
| <b>Entrez</b>            | <b>281249</b>                                                                                              |
| <b>Expiration Date</b>   | Please enquire.                                                                                            |