

Product Datasheet

Bovine Interleukin 17 (IL17) ELISA Kit (orb1207318)

Description	Bovine Interleukin 17 (IL17) ELISA Kit
Reactivity	Bovine
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	B3GDZ3
Sensitivity	11.2pg/mL
Entrez	282863
Expiration Date	Please enquire.