

Product Datasheet

Porcine Interleukin 15 (IL15) ELISA Kit (orb1207321)

Description	Porcine Interleukin 15 (IL15) ELISA Kit
Reactivity	Porcine
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q95253
Sensitivity	12.8pg/mL
Entrez	397683
Expiration Date	Please enquire.