

Product Datasheet

Bovine Interleukin 13 (IL13) ELISA Kit (orb1207324)

Description	Bovine Interleukin 13 (IL13) ELISA Kit
Reactivity	Bovine
Range	31.2-2000pg/ml
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9XSV9
Sensitivity	13.4pg/ml
Entrez	281247
Expiration Date	Please enquire.