

Product Datasheet

Porcine Interleukin 10 (IL10) ELISA Kit (orb1207327)

Description	Porcine Interleukin 10 (IL10) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q29055
Sensitivity	6.5pg/mL
Entrez	397106
Expiration Date	Please enquire.