

Product Datasheet

Porcine Insulin Like Growth Factor Binding Protein 1 (IGFBP1) ELISA Kit (orb1207336)

Description	Porcine Insulin Like Growth Factor Binding Protein 1 (IGFBP1) ELISA Kit
Reactivity	Porcine
Range	3.12-200ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	E0A8P1
Sensitivity	1.32ng/mL
Entrez	397270
Expiration Date	Please enquire.