

Product Datasheet

Porcine Interferon Gamma (IFN γ) ELISA Kit (orb1207342)

Description	Porcine Interferon Gamma (IFN γ) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P17803
Sensitivity	5.9pg/mL
Entrez	396991
Expiration Date	Please enquire.