

Product Datasheet

Porcine Fibroblast Growth Factor 2, Basic (FGF2) ELISA Kit (orb1207352)

Description	Porcine Fibroblast Growth Factor 2, Basic (FGF2) ELISA Kit
Reactivity	Porcine
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	F1RQN3
Sensitivity	41.78pg/mL
Entrez	397643
Expiration Date	Please enquire.