

## Product Datasheet

# Porcine Fibroblast Growth Factor 1, Acidic (FGF1) ELISA Kit (orb1207354)

|                          |                                                                                                      |
|--------------------------|------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Porcine Fibroblast Growth Factor 1, Acidic (FGF1) ELISA Kit                                          |
| <b>Reactivity</b>        | Porcine                                                                                              |
| <b>Range</b>             | 15.6-1000pg/mL                                                                                       |
| <b>Note</b>              | For research use only                                                                                |
| <b>Application notes</b> | sandwich                                                                                             |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Uniprot ID</b>        | <b>F1RM38</b>                                                                                        |
| <b>Sensitivity</b>       | 6.6pg/mL                                                                                             |
| <b>Entrez</b>            | <b>397497</b>                                                                                        |
| <b>Expiration Date</b>   | Please enquire.                                                                                      |