

Product Datasheet

Bovine Connective Tissue Growth Factor (CTGF) ELISA Kit (orb1207361)

Description	Bovine Connective Tissue Growth Factor (CTGF) ELISA Kit
Reactivity	Bovine
Range	62.5-4000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O18739
Sensitivity	30.23pg/mL
Entrez	281103
Expiration Date	Please enquire.