

Product Datasheet

Bovine Bone Morphogenetic Protein 2 (BMP2) ELISA Kit (orb1207365)

Description	Bovine Bone Morphogenetic Protein 2 (BMP2) ELISA Kit
Reactivity	Bovine
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	A5PJI9
Sensitivity	0.107ng/mL
Entrez	615037
Expiration Date	Please enquire.