

Product Datasheet

Mouse Plasminogen Activator, Urokinase Receptor (uPAR) ELISA Kit (orb1207388)

Description	Mouse Plasminogen Activator, Urokinase Receptor (uPAR) ELISA Kit
Reactivity	Mouse
Range	125-8000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P35456
Sensitivity	47pg/mL
Entrez	18793
Expiration Date	Please enquire.