

Product Datasheet

Rat Trefoil Factor 2 (TFF2) ELISA Kit (orb1207417)

Description	Rat Trefoil Factor 2 (TFF2) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q09030
Sensitivity	0.064ng/mL
Entrez	116592
Expiration Date	Please enquire.