

Product Datasheet

Mouse Superoxide Dismutase 3, Extracellular (SOD3) ELISA Kit (orb1207427)

Description	Mouse Superoxide Dismutase 3, Extracellular (SOD3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	009164
Sensitivity	0.064ng/mL
Entrez	20657
Expiration Date	Please enquire.