

Product Datasheet

Mouse S100A1 Calcium Binding Protein (S100A1) ELISA Kit (orb1207434)

Description	Mouse S100A1 Calcium Binding Protein (S100A1) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P56565
Sensitivity	28pg/mL
Entrez	20193
Expiration Date	Please enquire.