

Product Datasheet

Mouse P-Selectin Glycoprotein Ligand 1 (PSGL1) ELISA Kit (orb1207448)

Description	Mouse P-Selectin Glycoprotein Ligand 1 (PSGL1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q62170
Sensitivity	0.065ng/mL
Entrez	20345
Expiration Date	Please enquire.