

Product Datasheet

Mouse Paraoxonase 1 (PON1) ELISA Kit (orb1207455)

Description	Mouse Paraoxonase 1 (PON1) ELISA Kit
Reactivity	Mouse
Range	1.25-80ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P52430
Sensitivity	0.51ng/mL
Entrez	18979
Expiration Date	Please enquire.