

Product Datasheet

Mouse Phosphoenolpyruvate Carboxykinase 1, Soluble (PCK1) ELISA Kit (orb1207476)

Description	Mouse Phosphoenolpyruvate Carboxykinase 1, Soluble (PCK1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9Z2V4
Sensitivity	0.058ng/mL
Entrez	18534
Expiration Date	Please enquire.