

Product Datasheet

Mouse Osteocalcin (OC) ELISA Kit (orb1207487)

Description	Mouse Osteocalcin (OC) ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P86546
Sensitivity	23.58pg/mL
Entrez	12096
Expiration Date	Please enquire.