

Product Datasheet

Mouse Cross Linked N-Telopeptide Of Type I Collagen (NTXI) ELISA Kit (orb1207488)

Description	Mouse Cross Linked N-Telopeptide Of Type I Collagen (NTXI) ELISA Kit
Reactivity	Mouse
Range	3.12-200ng/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, urine or other biological fluids.
Sensitivity	0.96ng/mL
Expiration Date	Please enquire.