

Product Datasheet

Mouse Myosin Heavy Chain 2, Skeletal Muscle, Adult (MYH2) ELISA Kit (orb1207506)

Description	Mouse Myosin Heavy Chain 2, Skeletal Muscle, Adult (MYH2) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	G3UW82
Sensitivity	6.9pg/mL
Entrez	17882
Expiration Date	Please enquire.