

Product Datasheet

Mouse Mannose Associated Serine Protease 2 (MASP2) ELISA Kit (orb1207531)

Description	Mouse Mannose Associated Serine Protease 2 (MASP2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q91WP0
Sensitivity	0.067ng/mL
Entrez	17175
Expiration Date	Please enquire.