

Product Datasheet

Mouse Phospholipase A2, Lipoprotein Associated (LpPLA2) ELISA Kit (orb1207534)

Description	Mouse Phospholipase A2, Lipoprotein Associated (LpPLA2) ELISA Kit
Reactivity	Mouse
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q60963
Sensitivity	0.62ng/mL
Entrez	27226
Expiration Date	Please enquire.