

Product Datasheet

Mouse Heat Shock Protein 90kDa Beta 1 (HSP90b1) ELISA Kit (orb1207564)

Description	Mouse Heat Shock Protein 90kDa Beta 1 (HSP90b1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P08113
Sensitivity	0.064ng/mL
Entrez	22027
Expiration Date	Please enquire.