

Product Datasheet

Mouse Heme Oxygenase 1, Decycling (HO1) ELISA Kit (orb1207569)

Description	Mouse Heme Oxygenase 1, Decycling (HO1) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P14901
Sensitivity	6.3pg/mL
Entrez	15368
Expiration Date	Please enquire.