

Product Datasheet

Mouse Galectin 9 (GAL9) ELISA Kit (orb1207599)

Description	Mouse Galectin 9 (GAL9) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O08573
Sensitivity	29pg/mL
Entrez	16859
Expiration Date	Please enquire.